

## Pediatric Clubfoot Treatment



**P.E.A.care 2.0** Splint and  
**Clubycare** Boots

**Basko** Healthcare

Pieter Lieftinckweg 16 | 1505 HX Zaandam  
Tel.: +31 (0) 75 - 613 15 13 | Fax: +31 (0) 75 - 612 63 73  
E-mail: sales@basko.com | Internet: www.basko.com

[basko.com](http://basko.com)

# P.E.A.care 2.0 Static Quick Release

For the treatment of clubfoot with babies and toddlers Basko Healthcare provides orthotic components such as the P.E.A.care 2.0 splint and Clubycare boots. This combination in particular is used in the treatment of clubfoot according the Dr. Ponseti method.



## Features

- Quick Release mechanism provides easy donning and doffing of the boots
- Adjustable splint width
- Adjustable rotation angle setting (0° to 80°)
- Adjustable dorsalflexion (0° to 20°)
- Combinable with Clubycare boots
- Weight: 210 - 220 gram
- Material: durable XF Plastic

## Indications

- Arthrogryposis multiplex congenita (ICD10: Q74.3)
- Pes equinovarus congenitus (ICD10: Q66.0)

## Size chart - Static Quick Release

Art. no. "Static Quick Release"	Size	Splint width	Shoulder width
2850-1	1	195 - 275 mm	240 - 320 mm
2850-2	2	230 - 345 mm	280 - 390 mm

Measurement: The appropriate size is determined by the shoulder width.

# Clubycare Boots

For the treatment of clubfoot with babies and toddlers according the Dr. Ponseti method. Clubycare boots are made of soft leather and offer high comfort, tolerance and user compliance. The soles of the boots are easy to screw on to the foot plates of the P.E.A.care 2.0 splint. Clubycare boots are available with or without abduction in different sizes. The latter one comes either with lace closure or with velcro closure .

## Features

- Available with lace or velcro closure (905 - series)
- Padded and increased calcaneal cupping
- Open heel construction
- Open toe construction
- unit of measurement: **per piece**
- Elongated heel strap (stitched markers, unstretchable)
- Works perfectly in conjunction with P.E.A.care 2.0 splint
- Material: nappa leather; sole: rubber; lace: cotton; padding: EVA-foam

## Indications

- Arthrogryposis multiplex congenita (ICD10: Q74.3)
- Pes equinovarus congenitus (ICD10: Q66.0)



## Size chart - Clubycare Boots

Art. no. w/o abduction left with laces	Art. no. w/o abduction right with laces	Art. no. w/o abduction left with velcro	Art. no. w/o abduction right with velcro	Art. no. with abduction left with laces	Art. no. with abduction right with laces	Shoe size	Shoe length	Foot width
BA905L-12	BA905R-12	BA908L-12	BA908R-12	-	-	12	72 mm	40 mm
BA905L-13	BA905R-13	BA908L-13	BA908R-13	-	-	13	79 mm	42 mm
BA905L-14	BA905R-14	BA908L-14	BA908R-14	BA705L-14	BA705R-14	14	85 mm	43 mm
BA905L-15	BA905R-15	BA908L-15	BA908R-15	BA705L-15	BA705R-15	15	92 mm	45 mm
BA905L-16	BA905R-16	BA908L-16	BA908R-16	BA705L-16	BA705R-16	16	98 mm	46 mm
BA905L-17	BA905R-17	BA908L-17	BA908R-17	BA705L-17	BA705R-17	17	107 mm	47 mm
BA905L-18	BA905R-18	BA908L-18	BA908R-18	BA705L-18	BA705R-18	18	113 mm	48 mm
BA905L-19	BA905R-19	BA908L-19	BA908R-19	BA705L-19	BA705R-19	19	120 mm	50 mm
BA905L-20	BA905R-20	BA908L-20	BA908R-20	BA705L-20	BA705R-20	20	126 mm	51 mm
BA905L-21	BA905R-21	BA908L-21	BA908R-21	BA705L-21	BA705R-21	21	133 mm	52 mm
BA905L-22	BA905R-22	BA908L-22	BA908R-22	BA705L-22	BA705R-22	22	139 mm	54 mm
BA905L-23	BA905R-23	BA908L-23	BA908R-23	BA705L-23	BA705R-23	23	146 mm	55 mm
BA905L-24	BA905R-24	BA908L-24	BA908R-24	BA705L-24	BA705R-24	24	152 mm	56 mm
BA905L-25	BA905R-25	BA908L-25	BA908R-25	BA705L-25	BA705R-25	25	159 mm	58 mm
BA905L-26	BA905R-26	BA908L-26	BA908R-26	BA705L-26	BA705R-26	26	166 mm	59 mm
BA905L-27	BA905R-27	BA908L-27	BA908R-27	BA705L-27	BA705R-27	27	172 mm	60 mm
BA905L-28	BA905R-28	BA908L-28	BA908R-28	BA705L-28	BA705R-28	28	179 mm	62 mm
BA905L-29	BA905R-29	BA908L-29	BA908R-29	BA705L-29	BA705R-29	29	185 mm	63 mm
BA905L-30	BA905R-30	BA908L-30	BA908R-30	-	-	30	192 mm	64 mm

Sold per piece. An essential function of the Clubycare shoe is the heel strap, which ensures the heel is positioned correctly against the heel counter of the shoe. Therefore, sufficient growth reserve should be taken into account (max. 15 mm).